

Boston
Scientific



Quick Start Guide

Charging System

WAVEWRITER
ALPHA™

Spinal Cord Stimulator System



Ease of use is at the core of your SCS charging system

*The Charging System
aims to make life
as simple as
possible for you.*

An overview

Charger

Indicator Light



Power

+



Charging your device is as easy as 1–2–3!

- The charging system is **simple and easy to use**.
- The charging belt is **completely cordless**, so you can **enjoy full mobility** while charging your device.
- **You can decide** when and where you charge your device.



All you need to charge your rechargeable WaveWriter Alpha™ stimulator is:



A **power supply**



The **base station**



The **charger**



The **charging belt**

Charge your stimulator in 3 easy steps

1. Charge your charger

- Plug the base station into a power outlet and place the charger in the base station. A **yellow light** will appear on the charger. This means the **charging is in process**.



- Please note that as long as the **yellow light** shows, the charger is only **partially charged**. This means it can still be used to charge your stimulator but it may not be possible to recharge the stimulator fully.



- When a **green light** shows, the charger is **fully charged**. This means it is ready to charge your stimulator.




- If the indicator light is **not showing** while the charger is in the base station, there may be a **loose connection**. Check that the charger is securely placed in the base station.




2. Place the charger within your charging belt

- Once sufficiently charged, remove the charger from the base station. Before you place your charger in your charging belt, make sure to **switch your charger on**, using the circular **power button** in the middle of the charger.
- The **green indicator light** will show and the charger will **beep in a steady continuous rhythm**. Once the charger is correctly aligned with your stimulator, the **beeping will stop**.



-
- Place the **charger** into the pouch of your charger belt and place the belt around your waist.
 - Adjust the straps on the belt so they are loose enough for you to feel comfortable. Make sure the **mesh fabric of your belt is facing outwards**.



3. Position the belt to align the charger with your stimulator

- Move the belt slightly up and down and from side to side until the charger is **resting over your stimulator**. If you need to, you can adjust the length of the belt to the right size for you. When correctly positioned the **continuous single beeps will stop**.
- Check the charger occasionally to ensure it's **centered over your stimulator**. If not centered properly, charging times may vary.
- You know the stimulator is fully charged when you hear a **double beep for a total of 60 seconds** and see that the **green indicator light** has switched off. Although you don't have to charge your stimulator to 100% it can be helpful to know that your battery is fully charged.




Tip: Once your stimulator is charged, return your charger to the base station. This ensures it will be fully charged for the next time you need it. Then you're all set!

Making charging work for you

Top Tips

- **Choose a charging schedule that works for you:** You can decide whether to have shorter, more frequent charging sessions or longer, less frequent charging sessions. Pick a routine that best fits into your daily lifestyle and don't forget that you can easily charge while watching your favorite TV program or reading the morning paper. A standard charging session once a week generally takes **30–60 minutes**.
- **Check your battery status via your Freelink™ Remote Control:** As long as your charger isn't currently charging your stimulator, your Freelink™ remote will show your stimulator's battery status. Whenever this is low, your stimulator should be charged as soon as possible.
- **Keep your base station plugged in** so it is **always ready to charge** your charger when needed.



Learn more at RelieveMyPain.com.au

This document is not intended as a substitute for WaveWriter Alpha™ Spinal Cord Stimulator System Patient Manuals. The contents of the handbook must be carefully read and fully understood before using your charging system. If you have any questions, or need clarification on anything contained in this manual, please [contact your physician](#).

This material is for informational Purposes only and not meant for medical diagnosis. This information does not constitute medical or legal advice, and Boston Scientific makes no representation regarding the medical benefits included in this information. Boston Scientific strongly recommends that you consult with your physician on all matters pertaining to your health.

Talk to your healthcare professional about whether this product may be suitable for you as part of your overall plan to manage chronic pain.

The WaveWriter Alpha System is not a first-line treatment for chronic intractable pain.

Surgery is required in order to use the WaveWriter Alpha System and any surgical procedure carries risk.

Outcomes may vary for each patient.

Patients should consult their healthcare professional about factors that could impact their response.

THIS PRODUCT IS NOT AVAILABLE FOR PURCHASE BY THE GENERAL PUBLIC

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labelling supplied with each device or at www.IFU-BSCI.com. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries.

All trademarks are property of their respective owners.



Boston Scientific Pty Ltd
PO Box 332
Botany NSW 1455
Australia
Tel +61 2 8063 8100

© 2024 Boston Scientific Corporation or its affiliates. All rights reserved.

NM-1072104-AB NOV 2024